

1. T.M. Lecture Notes (Biology ID) (JND)

Ch. 1. Cell Biology

- ① Microscope
- ② Microbial Contamination
- ③ Heating Apparatus in Autoclave by Incubator
- ④ Centrifuge and Water Bath
- ⑤ Artificial Respiration

Course Lecturer: Dr. J. N. D.

✓ ① Microscope

It is a scientific instrument used to view an object that is too small to be seen by the naked eye.

Micro - small, scope - to look at or examine.
Scope - view hence the name microscope.

Types of Microscope

- a) Compound Microscope (Light Microscope)
- b) Electron Microscope
- c) Phase Contrast Microscope
- d) Interference Microscope
- e) Fluorescence Microscope
- f) Polarizing Microscope
- g) Dark field Microscope

The most commonly used is compound microscope. It is made of two converging lenses which are separated by distance which is known as eye piece. Magnification of compound microscope is objective lens x eye piece lens.

Parts of Microscope and Functions

- ① Eye piece: This is lens system near the eye which magnifies the primary image of the object so as to form a virtual image away from the eye point. With magnification of X10.

2. Rotatable head: To move the eye piece to the suitable position.

3. Arm: To carry the instrument comfortably.

4. Objective lenses: This is lens system which provide the image of the object. In general, the primary image they are form is inverted. Magnification of each are $\times 4$, $\times 10$, $\times 40$ and $\times 100$ at immersion.

5. Condenser: This is optical element which collect rays and collect to converge it through the stage through the small opening for the focus.

6. Base: This is the rectangular part which carries the illuminator, it is used for balancing on a smooth surface.

7. Illuminator: This is to provide light either from a sunlight or electricity and is housed at the base of the microscope.

8. Iris Diaphragm: This is mechanical device mounted within the condenser, it controlled the amount of light rays enter into objective lens parallel and to achieve optimum objective performance.

9. Nose Piece: This is the revolver that carries the objectives, it is hold to turn to the desired objective lens to be use.

10. Coarse Adjustment knob: This is knob using in moving the stage clockwise and anti-clockwise direction to achieve a focus. It is seen the side of the microscope.

11. Fine Adjustment knob: This is commonly found at the same place with coarse adjustment knob and is to give a fine and clearer focus.

CARE AND MAINTNANCE OF A MICROSCOPE

1. The mechanical and optical system must be handled with care.

2. Always carry the microscope with two hands.

(3)

1. Light of the arm and left hand under the base of the microscope.
2. Keep the microscope on a flat surface.
3. Wash all the microscope slide with soap and water.
4. After using (this should be neutral detergent) of various components.
5. Do not use organic solvents to wipe the surface of various components.
6. Never detach the eye-piece after use. Lower during adjustable knob, use to
7. No piece will be given while using objective.
8. Always keep the microscope inside the bag or box when using for safety.
9. Never disassemble the microscope for repair but reserve the task strictly for experienced technician.

MICROBIAL CONTAMINATION

Ever since the development of techniques to obtain micro organisms in pure culture, the susceptibility of such cultures to the unwanted growth of other microbes has been recognized.

This can be as a result of 'microbial contamination of culture media'.

In modern medicine and science, the importance of hand washing and the maintenance of a sterile operating theatre is taken for granted.

PREVENTION OF MICROBIOLOGICAL CONTAMINATION BEGINS IN THE LABORATORY

A variety of preventive procedures are common parts of an efficient microbiology laboratory. The use of sterile equipments and receptacles for liquid and solid growth media is a must. The prevention of contamination during the manipulations of microorganisms in the laboratory falls under the term 'aseptic technique'. Examples of aseptic technique include the disinfection

④
of work surfaces and the hands of the personnel before and after contact with microorganism and the flaming of the metal loops or rod used in transfer bacteria from one location to another.

First Aid measure for MICROBIC CONTAMINATION

In case of contact, wash contaminated skin thoroughly with antiseptic soap and water.

If wound contamination is suspected, seek immediate medical attention. In case of ingestion, seek immediate medical attention. In case of inhalation, seek immediate medical attention.

MICROTOME

Microtome is a machine used for cutting thin sections of tissue. This machine is used to study tissues with examination of tissue. It is made up of following components:

COMPONENTS OF MICROTOME

- ① A knife held in a firm support.
- ② A device for holding the piece of tissue to be cut.
- ③ A mechanism for advancing the piece of tissue across the knife blade or the knife across the piece.
- ④ A mechanism for advancing the piece of tissue across very small and accurate distance so that section of suitable thickness can be obtained.

TYPES OF MICROTOME

- 1 Cambridge Rocker Microtome.
- 2 Sledge Microtome with moving block.
- 3 Sledge Microtome with moving knife.
- 4 Rotary Microtome.
- 5 Freezing microtome.

(5)

STANDARD PROCEDURES FOR CARE OF MINOR INJURIES

1. General Instructions:

- A. Provide total care for those injured clearly within your capability to handle. All questionable cases should be referred to a local physician or to the emergency room.
- B. Wear disposable gloves and wash your hands thoroughly before and after treating an injury.

2. CARE OF INJURIES

A. Minor laceration (Cuts): Clean the wound with antiseptic solution and cover with appropriate dressing.

B. Severe laceration (Cuts): If free bleeding, apply pressure dressing and refer to doctor.

C. Minor cuts and abrasions: Apply antibiotic ointment and cover with appropriate dressing.

D. Minor abrasions (Scrapes): Gently cleanse area with appropriate antibiotic and cover with appropriate dressing.

E. Minor burns: Immerse burned area in cold water for 10-15 minutes, or until pain is relieved. Gently cleanse area, apply moist dressings and bandage loosely.

F. Foreign body in eye: If foreign body cannot be easily removed by flushing it with water, do not touch or try to remove it. Irrigate it with eye wash, refer to doctor.

CONTENTS OF THE FIRST AID BOX

- ✓ Name and phone number also address of the first aider or appointed person.
- ✓ Guidance leaflet.
- ✓ 20 individually wrapped sterile dressings.

10

dressings (assorted sizes) appropriate for the type of work

- ✓ 2 sterile eye pads with attached mounts
- ✓ 4 individually wrapped triaging wound bandages
- ✓ 6 safety pins
- ✓ 6 medium sized individually wrapped sterile UN-Medicated wound dressings (15cm x 12cm)
- ✓ 2 medium sized individually wrapped sterile UN-Medicated wound dressings (18cm x 18cm)
- ✓ One pair of disposable gloves
- ✓ Clinical waste bag with instruction for use

✓ LABORATORY SAFETY RULES

There are certain basic rules to be followed for safe working in the laboratories

- * Wear safety glasses and a fastened lab-coat at all times
- > Wear gloves if appropriate
- > Do not eat or drink or take food into the laboratory
- > Keep long hair tied back
- > Do not use a personal stereo or any other equipment involving headphones
- > Use all chemicals carefully and in a fume hood whenever possible
- > Dispose of solvents properly
- > Beware of glass especially if it is broken and then dispose of it safely
- > Keep your work area tidy and clean up any spills including water on the floor
- > If emergency alarm sound, make your work safe and leave the building as directed
- > Report any accident you have
- > Wash your hand thoroughly while leaving

The Laboratory

① Introduction To General Laboratory Techniques. (GLT-111)

Biology Option. NDI 2024/2025 Academic Session.

Lecture Note B1. Mr. M. CATH O'Riordan.

Introduction:

In a scientific investigation, there is an attitude that demands of one to explore his environment or to develop the relevant hypothesis in order to establish its acceptability or otherwise. This type of investigation can be carried out in a laboratory.

Laboratory:

A laboratory is a setting or place in order to try out specific information and test principles and theories in science.

It is thus, a place or room where scientific experiments are carried out. It is appropriately equipped with various items that are required to provide answers to interrogating questions in science.

Laboratory Hazards:

Working in a laboratory is characterized by certain specific features which in-turn make certain demands on the laboratory users.

There is need for special attention in order to adhere to necessary safety rules so as to minimize the occurrence of accident. There is also need to ensure that all experiments, reactions and related operations are carried out accurately.

Hazards:

Hazards can be defined as un-safe act that can bring about injury. Such injuries may be as a result of improper handling of glasswares, apparatus, reagents, chemicals and any other act of carelessness due to lack of obeying some elementary laboratory safety rules.

② Some SAFETY LABORATORY RULES.

There are certain basic rules to be followed for Safety Working in the laboratory. Such as :-

1. Do not eat or drink or take any food in the laboratory.
2. Remove if glasses especially if it is broken.
3. Keep your work area tidy and clean up any spills of water on the floor.
4. Keep long hair tied back and avoid using high shoes.
5. Be Careful and Concentrate on your work.
6. Wear your Lab-coat at all the time.
7. Wash your hands carefully before you leave the lab.
8. Report any accident you have to your supervisor.

CONTENTS OF FIRST AID BOX.

As a first Aider, you must have the following contents in your box.

First Aid is also defined as the first treatment given to the Patient before taking the Patient to the hospital or Doctor - if necessary.

1. Individual wrapped sterile adhesive dressing appropriate to the type of work.
2. Individual wrapped triangular bandages.
3. Safety Pins.
4. Unmedicated wound dressing.
5. Pairs of disposable gloves and nose mask.
6. Clinical Waste bag.
7. Packets of Rapid blade.
8. Water / cotton wool / iodine solution.
9. Sterilizer etc.

③ Microscope.

A Microscope is a scientific instrument used to magnify and observe some small objects that are too tiny to be seen with the naked eye.

It works by using lenses or electron beams to enlarge the image of the specimen. A microscope is widely used in Science, Medicine and Industry for research and diagnostic purposes.

There are two types of Microscope. Commonly used these are ① Light Compound Microscope ② Electron Microscope.

1. Light Compound Microscope

This type of Microscope uses visible light and lenses to magnify objects, commonly used in hospital and biology lab.

2. Electron Microscope

This type of Microscope uses electron beams for much higher magnification, allowing scientists to see fine structures like viruses and cell organelles.

Parts of Microscope.

A Microscope has several important parts which can be grouped into ① Optical ② Mechanical ③ Illuminating components.

1. Optical Components:

- Eye Piece:** - Is the lens you look through usually magnifies the image.
- Objective Lens:** - These are multiple lenses with different magnifications (eg. $\times 4$, $\times 10$, $\times 40$, $\times 100$ etc.).
- Revolving Nosepiece:** - It holds and rotates the objective lenses.

(4)

2. Mechanical Components:

- i Body tube: - This connects the eyepiece to the objective lenses
- ii Arm: - This supports the body tube and connects it to the base
- iii Stage: - Is the platform where the specimen slide is placed.
- iv Stage clip: - It holds the slide in place.
- v Coarse adjustment knob: - It moves the stage up and down for rough focusing.
- vi Fine adjustment knob: - It allows precise focusing or fine focusing.
- vii Base: - Is the bottom part that supports the microscope.

3. Illuminating Components:

- i Light source: - It provides light for viewing the specimen.
- ii Condenser lens: - It focuses light to the specimen.
- iii Diaphragm: - It controls the amount of light from passing through the specimen.

⑤

CARE AND MAINTENANCE OF A MICROSCOPE

Here are some important guidelines for care and maintenance of a microscope.

1. Handling and Storage:

- Always carry the microscope with two hands. One on the arm and one hand on the supporting base.
- Store the microscope in a clean, dry place, preferably covered when not in use.
- Keep the microscope on a stable vibration in free surface to prevent damage.

2.

Cleaning:

- Always use lens paper to clean the lenses, never use tissue or finger for cleaning.
- Clean the lenses with 70% alcohol using a cotton wool.
- Wipe the stage & body with a clean, dry cloth after use.

3.

Proper Use:

- Always start with the lowest objective lens eg $\times 4$ before switching to the higher objective lens.
- Use the coarse adjustment knob only with the low power objective lens to avoid crashing the lens in to the slide, and later on fine adjustment knob.
- Adjust the diaphragm and light source properly to get a clear image without straining your eyes.

4.

Preventive Maintenance:

- Regularly check and tighten any loose screws or knobs.
- Ensure that the light source works properly.
- Keep the microscope away from water and moisture.

Electron Microscope.

The electron microscope uses electrons to illuminate an object. Electrons have smaller wavelength than light, so they can resolve much smaller structures. Magnetic fields are used to create lenses that direct and focus the electrons. Since electrons are easily scattered by air molecules, the interior of the electron microscope must be used at a very high vacuum.

TYPES OF ELECTRON MICROSCOPE.

There are two types of electron microscope, these are: (1) Transmission electron microscope.

(2) Scanning electron microscope.

1. Transmission Electron Microscope: (TEM).

It is a powerful type of microscope that allows scientist to see ultra-thin slices of specimens in high details. Electron beam is directed on the object to be magnified. Some of the electrons are absorbed while others bounce off the specimen or pass through. It can magnify objects up to 1,00,000 times.

2. Scanning Electron Microscope: (SEM).

It focused beam of electrons to scan the surface of a specimen. The electron beam moves over the entire sample to create a magnified image of the surface of the object. It scans across electron beam on an image on to the screen of a television. It can also magnify object up to 100,000 times.

NOTE: Generally there are six (6) types of microscopes.

1. Compound Microscope
2. Electron microscope
3. Phase Contrast Microscope
4. Interference Microscope
5. Polarizing Microscope
6. Dark-field Microscope

Microtome.

Microtome is a Machine use for cutting thin structures or sections from pieces of tissue. This Machine is use in the study of histology which is the branch of medical science concerned with examination of tissue.

Many different varieties of Microtome are available, but all of them are made up of the following principal components.

1. A knife held in a firm support
2. A device for holding the piece of tissue firmly.
3. A mechanism for advancing the piece of tissue across the knife blade or the knife across the tissue.
4. A mechanism for advancing the piece of tissue across a very small and accurately measured distance so that section of suitable thickness can be obtained.

Types of Microtome.

There are five (5) types of microtome, these are,

1. Cambridge Rocker Microtome.
2. Sledge Microtome with moving knife.
3. Sledge Microtome with moving block.
4. Rotary Microtome.
5. Freezing Microtome.